

Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Device (Saliva)

Two ways to choose



20 Tests/kit

01 ▶



By Saliva collect cup

02 ▶



By Saliva collector

FEATURES:

- Specimen type: Saliva
- Testing time: 10 minutes
- Sensitivity: 93.94%
- Specificity: >99%


product name	Ref.	Size	Storage Temp
Novel Coronavirus(SARS-COV-2) Antigen Rapid Test Device (Saliva)	K590516D	20 Test/kit	2-30°C
Novel Coronavirus(SARS-COV-2) Antigen Rapid Test Device (Saliva)	K590516D	5 Test/kit	2-30°C
Novel Coronavirus(SARS-COV-2) Antigen Rapid Test Device (Saliva)	K590516D	1 Test/kit	2-30°C




The Novel Coronavirus (SARS-Cov-2)Antigen Rapid Test device(saliva)is an in vitro diagnostic test for the qualitative detection of novel coronavirus antigens in human saliva,using the rapid immunochromatographicmethod.The identification is based on the monoclonal antibodies specificfor the novel coronvirus antigen. It will provide information for clinicadoctors to prescribe correct medications.

DIRECTIONS FOR USE:


By Saliva collect cup



By Saliva collector



Test Results



HANGZHOU REALY TECH CO., LTD.

4th Floor,#12 Building,Eastern Medicine Town, Xiasha Economic&Technology Development,
310018 Hangzhou,Zhejiang,P.R.China

REALY 睿丽科技
Website: www.realytech.com